

## Product Datasheet

# Rat Lipolysaccharide Binding Protein ELISA Kit (orb1211128)

<b>Description</b>	Rat Lipolysaccharide Binding Protein, LBP ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	20-7000ng/l
<b>Tested Applications</b>	ELISA
<b>Target</b>	LBP
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	53,600
<b>Uniprot ID</b>	<a href="#">Q63313</a>
<b>Solubility (25°C)</b>	8000ng/L
<b>Sensitivity</b>	9.32ng/L
<b>Entrez</b>	<a href="#">29469</a>
<b>Expiration Date</b>	Please enquire.