

Product Datasheet

Human Vinculin ELISA Kit (orb1211136)

Description	Human Vinculin, VCL ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	VCL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	123,799
Uniprot ID	P18206
Solubility (25°C)	4000ng/L
Sensitivity	9.32ng/L
Entrez	7414
Expiration Date	Please enquire.