

Product Datasheet

Mouse Endostatin ELISA Kit (orb1211191)

Description	Mouse Endostatin, ES ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	ES
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	81,813
Uniprot ID	Q91V98
Solubility (25°C)	24ng/ml
Sensitivity	0.025ng/ml
Entrez	70445
Expiration Date	Please enquire.