

Product Datasheet

Mouse Proinsulin ELISA Kit (orb1211201)

Description	Mouse Proinsulin, PI ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	PI
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,458
Uniprot ID	P28798
Solubility (25°C)	480ng/L
Sensitivity	0.56ng/L
Entrez	14824
Expiration Date	Please enquire.