

Product Datasheet

Rat Macrophage Inflammatory Protein 1 Beta ELISA Kit (orb1211214)

Description	Rat Macrophage Inflammatory Protein 1 Beta, MIP-1B ELISA Kit
Reactivity	Rat
Range	1-400ng/L
Tested Applications	ELISA
Target	MIP-1beta
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,234
Uniprot ID	P50230
Solubility (25°C)	480ng/L
Sensitivity	0.54ng/L
Entrez	116637
Expiration Date	Please enquire.