

Product Datasheet

Canine Interleukin-2 Receptor Subunit Alpha ELISA Kit (orb1211226)

Description	Canine Interleukin-2 Receptor Subunit Alpha, IL2RA ELISA Kit
Reactivity	Canine
Range	30-9000ng/L
Tested Applications	ELISA
Target	IL2RA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,357
Uniprot ID	O62802
Solubility (25°C)	9600ng/L
Sensitivity	13.92ng/L
Entrez	403870
Expiration Date	Please enquire.