

## Product Datasheet

### Mouse Cardiac Troponin 1 ELISA Kit (orb1211229)

<b>Description</b>	Mouse Cardiac Troponin 1, CTN-1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	5-2000pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	cTn-1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	21,509
<b>Uniprot ID</b>	<a href="#">Q60753</a>
<b>Solubility (25°C)</b>	2400pg/ml
<b>Sensitivity</b>	2.47pg/ml
<b>Entrez</b>	<a href="#">13019</a>
<b>Expiration Date</b>	Please enquire.