

Product Datasheet

Rat Matrix Metalloproteinase 3 ELISA Kit (orb1211242)

Description	Rat Matrix Metalloproteinase 3, MMP-3 ELISA Kit
Reactivity	Rat
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	MMP-3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,428
Uniprot ID	P03957
Solubility (25°C)	64ng/ml
Sensitivity	0.07ng/ml
Entrez	171045
Expiration Date	Please enquire.