

Product Datasheet

Bovine Annexin 5 ELISA Kit (orb1211251)

Description	Bovine Annexin 5, ANX-5 ELISA Kit
Reactivity	Bovine
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	ANX-5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,089
Uniprot ID	P81287
Solubility (25°C)	40ng/ml
Sensitivity	0.022ng/ml
Entrez	281626
Expiration Date	Please enquire.