

Product Datasheet

Human Mucosal Addressin Cell Adhesion Molecule-1 ELISA Kit (orb1211253)

Description	Human Mucosal Addressin Cell Adhesion Molecule-1, MADCAM-1 ELISA Kit
Reactivity	Human
Range	20-3800ng/ml
Tested Applications	ELISA
Target	MAdCAM-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,155
Uniprot ID	Q13477
Solubility (25°C)	4000ng/ml
Sensitivity	9.19ng/ml
Entrez	8174
Expiration Date	Please enquire.