

Product Datasheet

Porcine Somatostatin ELISA Kit (orb1211270)

Description	Porcine Somatostatin, SS ELISA Kit
Reactivity	Porcine
Range	0.5-200ng/L
Tested Applications	ELISA
Target	SS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,429
Uniprot ID	P01248
Solubility (25°C)	240ng/L
Sensitivity	0.23ng/L
Entrez	396884
Expiration Date	Please enquire.