

## Product Datasheet

### Sheep Interleukin 8 ELISA Kit (orb1211363)

<b>Description</b>	Sheep Interleukin 8, IL-8 ELISA Kit
<b>Reactivity</b>	Sheep
<b>Range</b>	3-900ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL-8
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	11,292
<b>Uniprot ID</b>	<b>P36925</b>
<b>Solubility (25°C)</b>	960ng/L
<b>Sensitivity</b>	1.62ng/L
<b>Entrez</b>	<b>443418</b>
<b>Expiration Date</b>	Please enquire.