

Product Datasheet

Rat Neuropeptide Y Receptor Type 5 ELISA Kit (orb1211419)

Description	Rat Neuropeptide Y Receptor Type 5, NPY5R ELISA Kit
Reactivity	Rat
Range	20-6000ng/L
Tested Applications	ELISA
Target	NPY5R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,408
Uniprot ID	Q63634
Solubility (25°C)	6400ng/L
Sensitivity	9.65ng/L
Entrez	25340
Expiration Date	Please enquire.