

Product Datasheet

Rat Metallothionein-2 ELISA Kit (orb1211434)

Description	Rat Metallothionein-2, MT-2 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	MT2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	6,145
Uniprot ID	P04355
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	689415
Expiration Date	Please enquire.