

Product Datasheet

Rat Matrix Metalloproteinase-16 ELISA Kit (orb1211437)

Description	Rat Matrix Metalloproteinase-16, MMP16 ELISA Kit
Reactivity	Rat
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	MMP16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,624
Uniprot ID	035548
Solubility (25°C)	96ng/mL
Sensitivity	0.141ng/mL
Entrez	65205
Expiration Date	Please enquire.