

Product Datasheet

Rat Lactadherin ELISA Kit (orb1211441)

Description	Rat Lactadherin, MFG8 ELISA Kit
Reactivity	Rat
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	MFG8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,413
Uniprot ID	P70490
Solubility (25°C)	240ng/mL
Sensitivity	0.242ng/mL
Entrez	25277
Expiration Date	Please enquire.