

Product Datasheet

Rat Fibroblast Growth Factor 1 ELISA Kit (orb1211452)

Description	Rat Fibroblast Growth Factor 1, FGF1 ELISA Kit
Reactivity	Rat
Range	3-900ng/L
Tested Applications	ELISA
Target	FGF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,418
Uniprot ID	P61149
Solubility (25°C)	960ng/L
Sensitivity	1.79ng/L
Entrez	25317
Expiration Date	Please enquire.