

Product Datasheet

Rat Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha ELISA Kit (orb1211457)

Description	Rat Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha, PCG-1A ELISA Kit
Reactivity	Rat
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	PPARGC1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,622
Uniprot ID	Q9QYK2
Solubility (25°C)	32ng/ml
Sensitivity	0.023ng/ml
Entrez	83516
Expiration Date	Please enquire.