

## Product Datasheet

### Rat Glutaredoxin-1 ELISA Kit (orb1211468)

<b>Description</b>	Rat Glutaredoxin-1, GLRX ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.5-200ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	GLRX
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	11,879
<b>Uniprot ID</b>	<a href="#">Q9ESH6</a>
<b>Solubility (25°C)</b>	240ng/mL
<b>Sensitivity</b>	0.28ng/mL
<b>Entrez</b>	<a href="#">64045</a>
<b>Expiration Date</b>	Please enquire.