

Product Datasheet

Rat Galanin Receptor Type 1 ELISA Kit (orb1211477)

Description	Rat Galanin Receptor Type 1, GALR1 ELISA Kit
Reactivity	Rat
Range	20-4500ng/L
Tested Applications	ELISA
Target	GALR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,862
Uniprot ID	Q62805
Solubility (25°C)	4800ng/L
Sensitivity	9.03ng/L
Entrez	50577
Expiration Date	Please enquire.