

Product Datasheet

Rat Fatty Acid Synthase ELISA Kit (orb1211484)

Description	Rat Fatty Acid Synthase, FAS ELISA Kit
Reactivity	Rat
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	FASN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	272,650
Uniprot ID	P12785
Solubility (25°C)	240ng/mL
Sensitivity	0.27ng/mL
Entrez	50671
Expiration Date	Please enquire.