

Product Datasheet

Rat Cx3C Chemokine Receptor 1 ELISA Kit (orb1211495)

Description	Rat Cx3C Chemokine Receptor 1, CX3CR1 ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	CX3CR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,327
Uniprot ID	P35411
Solubility (25°C)	3200ng/L
Sensitivity	7.62ng/L
Entrez	171056
Expiration Date	Please enquire.