

Product Datasheet

Rat Decorin ELISA Kit (orb1211539)

Description	Rat Decorin, DCN ELISA Kit
Reactivity	Rat
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	DCN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,805
Uniprot ID	Q01129
Solubility (25°C)	32ng/ml
Sensitivity	0.036ng/ml
Entrez	29139
Expiration Date	Please enquire.