

Product Datasheet

Rat Adisintegrin and Metalloproteinase with Thrombospondin Motifs 4 ELISA Kit (orb1211562)

Description	Rat Adisintegrin and Metalloproteinase with Thrombospondin Motifs 4, ADAMTS4 ELISA Kit
Reactivity	Rat
Range	5-300ng/ml
Tested Applications	ELISA
Target	ADAMTS4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,385
Uniprot ID	Q9ESP7
Solubility (25°C)	320ng/ml
Sensitivity	2.02ng/ml
Entrez	66015
Expiration Date	Please enquire.