

Product Datasheet

Rat Peroxiredoxin 5 ELISA Kit (orb1211567)

Description	Rat Peroxiredoxin 5, PRDX5 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	PRDX5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,179
Uniprot ID	Q9R063
Solubility (25°C)	24ng/ml
Sensitivity	0.022ng/ml
Entrez	113898
Expiration Date	Please enquire.