

Product Datasheet

Rat Transferrin Receptor Protein 1 ELISA Kit (orb1211576)

Description	Rat Transferrin Receptor Protein 1, TFR1 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	TFRC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	85,876
Uniprot ID	Q99376
Solubility (25°C)	24ng/ml
Sensitivity	0.021ng/ml
Entrez	64678
Expiration Date	Please enquire.