

Product Datasheet

Rat Lactoylglutathione Lyase ELISA Kit (orb1211580)

| | |
|----------------------------|---|
| Description | Rat Lactoylglutathione Lyase, GLO1 ELISA Kit |
| Reactivity | Rat |
| Range | 0.5-200ng/ml |
| Tested Applications | ELISA |
| Target | GLO1 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 20,820 |
| Uniprot ID | Q6P7Q4 |
| Solubility (25°C) | 240ng/ml |
| Sensitivity | 0.29ng/ml |
| Entrez | 294320 |
| Expiration Date | Please enquire. |