

## Product Datasheet

# Rat Matrix Metalloproteinase 14 ELISA Kit (orb1211584)

<b>Description</b>	Rat Matrix Metalloproteinase 14, MMP14 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	30-7000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	MMP14
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	66,080
<b>Uniprot ID</b>	<a href="#">Q10739</a>
<b>Solubility (25°C)</b>	8000ng/L
<b>Sensitivity</b>	14.96ng/L
<b>Entrez</b>	<a href="#">81707</a>
<b>Expiration Date</b>	Please enquire.