

Product Datasheet

Rat Insulin-Degrading Enzyme ELISA Kit (orb1211623)

Description	Rat Insulin-Degrading Enzyme, IDE ELISA Kit
Reactivity	Rat
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	IDE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	117,710
Uniprot ID	P35559
Solubility (25°C)	80ng/ml
Sensitivity	0.12ng/ml
Entrez	25700
Expiration Date	Please enquire.