

Product Datasheet

Rat Albumin ELISA Kit (orb1211647)

Description	Rat Albumin ELISA Kit
Reactivity	Rat
Range	1-300ug/ml
Tested Applications	ELISA
Target	ALB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,731
Uniprot ID	P02770
Solubility (25°C)	320ug/ml
Sensitivity	0.47ug/ml
Entrez	24186
Expiration Date	Please enquire.