

## Product Datasheet

### Rat Afamin ELISA Kit (orb1211652)

<b>Description</b>	Rat Afamin, AFM ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	1-300ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	AFM
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	69,335
<b>Uniprot ID</b>	<b>P36953</b>
<b>Solubility (25°C)</b>	320ng/ml
<b>Sensitivity</b>	0.42ng/ml
<b>Entrez</b>	<b>282708</b>
<b>Expiration Date</b>	Please enquire.