

Product Datasheet

Rat Gamma Glutamyl Transpeptidase 1 ELISA Kit (orb1211667)

Description	Rat Gamma Glutamyl Transpeptidase 1, GGT1 ELISA Kit
Reactivity	Rat
Range	2-600mIU/mL
Tested Applications	ELISA
Target	GGT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,610
Uniprot ID	P07314
Solubility (25°C)	640mIU/mL
Sensitivity	1.13mIU/mL
Entrez	116568
Expiration Date	Please enquire.