

Product Datasheet

Rat Secreted Frizzled Related Protein 4 ELISA Kit (orb1211677)

Description	Rat Secreted Frizzled Related Protein 4, SFRP4 ELISA Kit
Reactivity	Rat
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	SFRP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,762
Uniprot ID	Q9JLS4
Solubility (25°C)	12.8ng/ml
Sensitivity	0.025ng/ml
Entrez	89803
Expiration Date	Please enquire.