

## Product Datasheet

# Rat Cxc-Chemokine Receptor 4 ELISA Kit (orb1211696)

<b>Description</b>	Rat Cxc-Chemokine Receptor 4, CXCR4 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.1-30ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	CXCL4
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	11,287
<b>Uniprot ID</b>	<b>P06765</b>
<b>Solubility (25°C)</b>	32ng/mL
<b>Sensitivity</b>	0.06ng/mL
<b>Entrez</b>	<b>360918</b>
<b>Expiration Date</b>	Please enquire.