

Product Datasheet

Rat Oxidized Low-density Lipoprotein Receptor 1 ELISA Kit (orb1211697)

Description	Rat Oxidized Low-density Lipoprotein Receptor 1, LOX-1 ELISA Kit
Reactivity	Rat
Range	0.5-60ng/ml
Tested Applications	ELISA
Target	OLR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,890
Uniprot ID	O70156
Solubility (25°C)	64ng/ml
Sensitivity	0.32ng/ml
Entrez	140914
Expiration Date	Please enquire.