

Product Datasheet

Rat Transferrin Receptor Protein 2 ELISA Kit (orb1211703)

| | |
|----------------------------|--|
| Description | Rat Transferrin Receptor Protein 2, TFR2 ELISA Kit |
| Reactivity | Rat |
| Range | 5-1000ug/ml |
| Tested Applications | ELISA |
| Target | TFR2 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 88,328 |
| Uniprot ID | B2GUY2 |
| Solubility (25°C) | 1280µg/ml |
| Sensitivity | 2.42µg/ml |
| Entrez | 288562 |
| Expiration Date | Please enquire. |