

Product Datasheet

Rat Glucagon-like Peptide 1 Receptor ELISA Kit (orb1211718)

Description	Rat Glucagon-like Peptide 1 Receptor, GLP1R ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	GLP1R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,877
Uniprot ID	P32301
Solubility (25°C)	3200ng/L
Sensitivity	7.11ng/L
Entrez	25051
Expiration Date	Please enquire.