

Product Datasheet

Rat Sorbitol Dehydrogenase ELISA Kit (orb1211720)

Description	Rat Sorbitol Dehydrogenase, SDH ELISA Kit
Reactivity	Rat
Range	0.1-40mIU/mL
Tested Applications	ELISA
Target	SORD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,235
Uniprot ID	P27867
Solubility (25°C)	48mIU/mL
Sensitivity	0.056mIU/mL
Entrez	24788
Expiration Date	Please enquire.