

Product Datasheet

Rat Kisspeptin 1 ELISA Kit (orb1211723)

Description	Rat Kisspeptin 1, KISS-1 ELISA Kit
Reactivity	Rat
Range	0.06-18ng/mL
Tested Applications	ELISA
Target	KISS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,188
Uniprot ID	Q7TSB7
Solubility (25°C)	19.2ng/mL
Sensitivity	0.034ng/mL
Entrez	289023
Expiration Date	Please enquire.