

Product Datasheet

Rat Sclerostin Domain-containing Protein 1 ELISA Kit (orb1211739)

Description	Rat Sclerostin Domain-containing Protein 1, SOSTDC1 ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	SOSTDC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,136
Uniprot ID	Q642G2
Solubility (25°C)	3200ng/L
Sensitivity	7.37ng/L
Entrez	266803
Expiration Date	Please enquire.