

## Product Datasheet

### Rat Aquaporin 2 ELISA Kit (orb1211783)

<b>Description</b>	Rat Aquaporin 2, AQP-2 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	3-900ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	AQP2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	28,931
<b>Uniprot ID</b>	<a href="#">P34080</a>
<b>Solubility (25°C)</b>	960ng/mL
<b>Sensitivity</b>	1.56ng/mL
<b>Entrez</b>	<a href="#">25386</a>
<b>Expiration Date</b>	Please enquire.