

Product Datasheet

Rat Monocyte Chemotactic Protein 2 ELISA Kit (orb1211846)

Description	Rat Monocyte Chemotactic Protein 2, MCP-2 ELISA KI
Reactivity	Rat
Range	2-600ng/L
Tested Applications	ELISA
Target	CCL2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,460
Uniprot ID	P14844
Solubility (25°C)	800ng/L
Sensitivity	1.01ng/L
Entrez	24770
Expiration Date	Please enquire.