

Product Datasheet

Rat Interleukin 6 ELISA Kit (orb1211853)

Description	Rat Interleukin 6, IL-6 ELISA Kit
Reactivity	Rat
Range	0.1-40ng/L
Tested Applications	ELISA
Target	IL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,357
Uniprot ID	P20607
Solubility (25°C)	48ng/L
Sensitivity	0.052ng/L
Entrez	24498
Expiration Date	Please enquire.