

Product Datasheet

Rat Interleukin 33 ELISA Kit (orb1211855)

Description	Rat Interleukin 33, IL-33 ELISA Kit
Reactivity	Rat
Range	1-400ng/L
Tested Applications	ELISA
Target	IL33
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,486
Uniprot ID	Q66H70
Solubility (25°C)	480ng/L
Sensitivity	0.56ng/L
Entrez	361749
Expiration Date	Please enquire.