

Product Datasheet

Rat Orexin A ELISA Kit (orb1211863)

Description	Rat Orexin A ELISA Kit
Reactivity	Rat
Range	7-1500ng/L
Tested Applications	ELISA
Target	HCRT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	13,645
Uniprot ID	O55232
Solubility (25°C)	1600ng/L
Sensitivity	3.15ng/L
Entrez	25723
Expiration Date	Please enquire.