

Product Datasheet

Rat Toll-like Receptor 4 ELISA Kit (orb1211867)

Description	Rat Toll-like Receptor 4, TLR4 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	TLR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,072
Uniprot ID	Q9QX05
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	29260
Expiration Date	Please enquire.