

Product Datasheet

Rabbit Apolipoprotein A1 ELISA Kit (orb1211883)

| | |
|----------------------------|---|
| Description | Rabbit Apolipoprotein A1, APO-A1 ELISA Kit |
| Reactivity | Rabbit |
| Range | 2-700ug/ml |
| Tested Applications | ELISA |
| Target | APOA1 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 30,591 |
| Uniprot ID | P09809 |
| Solubility (25°C) | 800µg/ml |
| Sensitivity | 1.21µg/ml |
| Entrez | 100009253 |
| Expiration Date | Please enquire. |