

Product Datasheet

Rabbit Insulin-like Growth Factors 1 ELISA Kit (orb1211884)

| | |
|----------------------------|---|
| Description | Rabbit Insulin-like Growth Factors 1, IGF-1 ELISA Kit |
| Reactivity | Rabbit |
| Range | 5-1000ng/ml |
| Tested Applications | ELISA |
| Target | IGF1 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 16,092 |
| Uniprot ID | Q95222 |
| Solubility (25°C) | 1280ng/ml |
| Sensitivity | 2.25ng/ml |
| Entrez | 100008668 |
| Expiration Date | Please enquire. |