

Product Datasheet

Rabbit Endothelin 1 ELISA Kit (orb1211891)

| | |
|----------------------------|--|
| Description | Rabbit Endothelin 1, ET-1 ELISA Kit |
| Reactivity | Rabbit |
| Range | 3-900ng/L |
| Tested Applications | ELISA |
| Target | EDN1 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 22,828 |
| Uniprot ID | P29560 |
| Solubility (25°C) | 960ng/L |
| Sensitivity | 1.63ng/L |
| Entrez | 100009270 |
| Expiration Date | Please enquire. |