

Product Datasheet

Porcine Interferon Alpha ELISA Kit (orb1211913)

Description	Porcine Interferon Alpha, Beta Receptor 1, IFNAR1 ELISA Kit
Reactivity	Porcine
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	IFNAR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,221
Uniprot ID	Q764M8
Solubility (25°C)	48ng/mL
Sensitivity	0.055ng/mL
Entrez	396658
Expiration Date	Please enquire.